

LED Filter Set Selection Guide



Filter Set

More filters are selectable as your experiment needs

<i>Product No.</i>	<i>Model No.</i>	<i>Filter Diameter</i>	<i>Filter(EX)</i>	<i>Dichroic Mirror</i>	<i>Filter(EM)</i>
49028	Custom DAPI	25mm	AT375/28x	T425lpxr	ET460/50m
49000	DAPI	25mm	AT350/50x	T4000lp	ET460/50m
49001	EYFP/Cerulean	25mm	AT435/20x	AT455DC	AT480/30m
49002	EGFP (FITC/Cy2)	25mm	ET470/40x	T495lpxr	ET525/50m
49003	EYFP	25mm	ET500/20x	T515lp	ET535/30m
49004	CY3/TRITC	25mm	ET545/25x	T565lpxr	ET605/70m
49005	DSRed (TRITC/Cy3)	25mm	ET545/30x	T570lp	ET620/60m
49008	mCherry, Texas Red	25mm	ET560/40x	T585lpxr	ET630/75m
49009	Cy5 Narrow Excitation	25mm	ET640/30x	T660lpxr	ET690/50m
49006	Cy5	25mm	ET620/60x	T660lpxr	ET700/75m